


STAGEFLEXER® MEMBRANES

Flexible growth surface for the StageFlexer and FlexFlow devices

- Individually packaged 43 mm silicone rubber membranes.
- Membranes come in sterile culture dishes.
- Clear for direct viewing of cells optically.
- Viscoelastic under mechanical loading.
- Same growth surfaces as the BioFlex culture plates.
- Covalently bonded surfaces: Amino, Collagen (Type I or IV) Elastin, ProNectin (RGD), Laminin (YIGSR).



Ordering Information	
<i>Catalog Number</i>	<i>Product/Item</i>
SFM-U	StageFlexer Membrane - Untreated
SFM-A	StageFlexer Membrane - Amino
SFM-C	StageFlexer Membrane - Collagen Type I
SFM-C/IV	StageFlexer Membrane - Collagen Type IV
SFM-E	StageFlexer Membrane - Elastin
SFM-P	StageFlexer Membrane - ProNectin
SFM-L	StageFlexer Membrane - Laminin