



FLEXCELL® FLEXLINK®

Tension or Compression

(cat. nos. FX-5000TFL or FX-5000CFL)

Run multiple regimens at one time!



- Strain unit module that links with a FX-5000 system computer allowing the user to simultaneously operate two different regimens on two different baseplates.
- Connects in a series to operate up to four FlexLinks from a central computer.
- FlexLink may require an additional vacuum or positive pressure source depending on the intensity of the regimens.



Figure 10. Image of FlexLink Tension System controller.

FX-5000 FlexLink includes:

- Air Handling: Valves and Pressure Transducers
- Accessory package:
 - One Baseplate and four gaskets
 - Tubing and quick disconnects