

# SHEL LAB *BacBASiC*

Environmental Work Station & Incubator  
For Anaerobic Applications



Ergonomic, "Gun  
Drop"\* Arm Port  
Doors offer  
maximum reach  
& comfort!

\*ARM PORT  
DOOR DESIGN  
PATENT  
PENDING



Removable  
viewing shield for  
introduction of  
equipment,  
cleaning,  
decontamination  
or maintenance.

The only thing "basic" about the SHEL LAB *BacBASiC* is its very reasonable price.

The *BacBASiC* is the ideal, entry level chamber for performance of truly anaerobic microbiology. Move beyond jars and the unnecessary exposure of your anaerobes to deadly oxygen.

*Already have a chamber?* Then the *BacBASiC* offers the perfect auxiliary chamber option for over-flow work, extra plate storage and emergency back-up for your primary chamber!

Constant thinking.

**SL**  
SHEL LAB

A Product of  
SHELDON MFG. INC.

# SHEL LAB BacBASIC

The SHEL LAB **BacBASIC** offers:

- An ample, ergonomic, user-friendly work station with a minimal, space-saving exterior footprint.
- Excellent lighting and visibility.
- A quick-entry, gloveless sleeve system.
- A rapid\* Auto-Purge Pass Box for introducing Petri plates, specimens and supplies. (\*Purge time less than 50 seconds.)
- Incubated storage capacity for up to 150 Petri plates.
- Incubation temperature range from ambient temperature to 42°C.
- Easy-Change anatox and catalyst packs and the catalyst can be used for up to six months before replacement or reactivation.



*Simple Push  
Button Operation*

*Easy Change Anatox  
& Catalyst Packs*

*Central Pass  
Through Box*

*Interior Auto-  
Purge Button*

For more information on the  
**SHEL LAB BacBASIC:**  
Call: (503) 640-3000  
Go To: [www.shellab.com](http://www.shellab.com)

## SHEL LAB **BacBASIC**

Exterior Dimensions (w x d x h)	34.75 in x 29.5 in x 29.25 in 88.3 cm x 74.9 cm x 74.3 cm
Interior Dimensions (w x d x h)	32.5 in x 19.25 in x 22 in 82.6 cm x 48.6 cm x 55.9 cm
Weight	230 lbs / 104.3 kg
Temperature Range	Ambient to 42°C
Temperature Uniformity	+/- 0.5°C

**Dunn**  
**Labortechnik GmbH**



Thelenberg 6 · 53567 Asbach · Tel. 0 26 83 / 4 30 94 · Fax 0 26 83 / 4 27 76  
E-mail: [info@dunnlab.de](mailto:info@dunnlab.de) · Internet: [www.dunnlab.de](http://www.dunnlab.de)

Constant thinking.

**SL**  
**SHEL LAB**  
A Product of  
SHELDON MFG. INC.