

Erlenmeyer-, Fernbach- and Shake Flasks from Borosilicate Glass

- complete price list on request -

Erlenmeyer Flasks (graduated)

For general laboratory use such as media preparation, storage and cultivation of microorganisms

Narrow Mouth with Break-Resistant Lip

Cat. No.	Capacity (ml)	Stopper Size	Ø x Height (mm)	Cs.
2504-00025	25	0	39 x 66	24
2504-00050	50	1	50 x 80	24
2504-00125	125	5	66 x 115	12
2504-00250	250	6	80 x 140	12
2504-00300	300	6	88 x 148	12
2504-00500	500	7	100 x 183	12
2504-01000	1000	9	128 x 223	6
2504-02000	2000	10	158 x 290	3
2504-03000	3000	10	184 x 305	3
2504-04000	4000	10	205 x 355	3
2504-06000	6000	10	237 x 395	3



Narrow Mouth

Cat. No.	Capacity (ml)	Stopper Size	Ø x Height (mm)	Cs.
2505-00125	125	5	66 x 107	12
2505-00250	250	6	80 x 133	12
2505-00500	500	7	100 x 173	6
2505-01000	1000	9	128 x 214	6



Narrow Mouth with Phenolic Rubber Lined Screw Cap

- Long-skirted screw caps protect against contamination

Cat. No.	Capacity (ml)	Cap (mm)	Ø x Height (mm)	Cs.
2514-00025	25	18-415	40 x 65	24
2514-00050	50	24-415	50 x 85	24
2514-00125	125	24-415	66 x 125	12
2514-00250	250	38-415	80 x 143	12
2514-00300	300	38-415	88 x 150	12
2514-00500	500	38-415	100 x 190	12
2514-01000	1000	38-415	128 x 220	6
2514-02000	2000	38-415	158 x 275	3
2514-03000	3000	38-415	184 x 320	3
2514-04000	4000	38-415	205 x 350	3
2514-06000	6000	38-415	235 x 390	3



Wide Mouth with GL 45 Screw Cap and Pouring Ring

- With 45 mm Polypropylene Linerless Screw Cap and Pouring Ring

Cat. No.	Capacity (ml)	Cap (mm)	Ø x Height (mm)	Cs.
2515-00250	250	45	80 x 117	12
2515-00500	500	45	100 x 156	12
2515-01000	1000	45	128 x 200	6
2515-02000	2000	45	158 x 260	3
2515-04000	4000	45	205 x 325	1
2515-06000	6000	45	235 x 370	1



Fernbach Flasks

For culturing organism requiring a large surface to volume ratio

DeLong® Style Neck, Ø 38 mm

- Without Baffles

Cat. No.	Capacity (ml)	Ø x Height (mm)	Cs.
2550-01800	1800	200 x 175	3
2550-02800	2800	205 x 250	3



Wide Mouth with Break-Resistant Lip

- With three Side Baffles

Cat. No.	Capacity (ml) and Description	Ø x Height (mm)	Cs.
2551-01800	1800, 3 Side Baffles	200 x 157	3
2551-02800	2800, 3 Side Baffles	210 x 225	3



DeLong® Style Neck, Ø 38 mm

- With three Bottom Baffles, or three Bottom and three Side Baffles

Cat. No.	Capacity (ml) and Description	Ø x Height (mm)	Cs.
2554-01800	1800, 3 Bottom Baffles	200 x 175	3
2554-02800	2800, 3 Bottom Baffles	210 x 250	3
2554-02833	2800, 6 Baffles	200 x 250	3



Wide Mouth with GL 45 Screw Cap and Pouring Ring

- 45 mm Polypropylene Screw Cap
- With three Bottom Baffles, or three Bottom and three Side Baffles

Cat. No.	Capacity (ml) and Description	Ø x Height (mm)	Cs.
2552-01800	1800, 3 Bottom Baffles	200 x 157	1
2552-02800	2800, 3 Bottom Baffles	210 x 225	1
2552-02833	2800, 6 Baffles	210 x 225	1



Wide Mouth with GL 45 Membrane Screw Cap and Pouring Ring

- 45 mm Polypropylene Screw Cap with 0,2 µm PTFE Membrane
- With three Bottom Baffles, or three Bottom and three Side Baffles

Cat. No.	Capacity (ml) and Description	Ø x Height (mm)	Cs.
2552-18000	1800, no Baffles	200 x 157	1
2552-28000	2800, no Baffles	210 x 225	1
2552-01803	1800, 3 Bottom Baffles	200 x 157	1
2552-02803	2800, 3 Bottom Baffles	210 x 225	1
2552-62180	1800 ml, 6 Baffles	200 x 157	1
2552-62280	1800 ml, 6 Baffles	210 x 225	1
5637-00046	GL 45, Polypropylene, 0,2 µm PTFE Membrane Cap		5

Shake Flasks

For preliminary fermentation and bacterial growth studies

DeLong® Style Neck, three Side Baffles

Cat. No.	Capacity (ml)	Closures (mm)	Ø x Height (mm)	Cs.
2540-00050	50	18	50 x 105	12
2540-00125	125	25	66 x 136	12
2540-00250	250	38	80 x 160	12
2540-00300	300	38	88 x 170	12
2540-00500	500	38	100 x 205	12
2540-01000	1000	38	128 x 235	6
2540-02000	2000	38	158 x 290	3



Break-Resistant Lip, three Side Baffles

Cat. No.	Capacity (ml)	Ø x Height (mm)	Cs.
2541-00125	125	66 x 116	12
2541-00250	250	80 x 140	12
2541-01000	1000	128 x 220	6
2541-02000	2000	158 x 290	3



Break-Resistant Lip, three Extra Deep Side Baffles

Cat. No.	Capacity (ml)	Ø x Height (mm)	Cs.
2542-00250	250	80 x 140	12
2542-00300	300	88 x 150	12
2542-00500	500	100 x 183	12
2542-01000	1000	128 x 220	6
2542-02000	2000	158 x 290	3



DeLong® Style Neck, three Extra Deep Side Baffles


Cat. No.	Capacity (ml)	Closures (mm)	Ø x Height (mm)	Cs.
2543-00050	50	18	50 x 105	12
2543-00125	125	25	66 x 136	12
2543-00250	250	38	80 x 160	12
2543-00300	300	38	88 x 170	12
2543-00500	500	38	100 x 205	12
2543-01000	1000	38	128 x 235	6
2543-02000	2000	38	158 x 295	3
2543-03000	3000	38	184 x 340	3
2543-04000	4000	38	205 x 370	3



DeLong® Style Neck, six Baffles

- With three Bottom Baffles, or three Bottom and three Side Baffles

Cat. No.	Capacity (ml)	Closures (mm)	Ø x Height (mm)	Cs.
2543-60250	250	38	80 x 160	12
2543-60500	500	38	100 x 205	12
2543-61000	1000	38	128 x 235	6
2543-62000	2000	38	158 x 295	3
2543-63000	3000	38	184 x 340	3
2543-64000	4000	38	205 x 370	3



Wide Mouth with GL 45 PP Membrane Screw Cap and Pouring Ring

- 45 mm Polypropylene Screw Cap with 0,2 µm PTFE Membrane (5637-00046)
- With three Bottom Baffles, or three Bottom and three Side Baffles

Cat. No.	Capacity (ml)	Closures (mm)	Ø x Height (mm)	Cs.
2545-00250	250	45	80 x 117	12
2545-00500	500	45	100 x 156	12
2545-01000	1000	45	128 x 200	6
2545-02000	2000	45	158 x 260	3



DeLong® Style Neck, four Extra Deep Bottom Baffles, with Constriction

Cat. No.	Capacity (ml)	Closures (mm)	Ø x Height (mm)	Cs.
2547-00250	250	38	81 x 160	12
2547-00500	500	38	100 x 203	12
2547-01000	1000	38	128 x 235	6



DeLong® Style Neck with Stainless Steel Closures, without Baffles

Cat. No.	Capacity (ml)	Closures (mm)	Ø x Height (mm)	Cs.
2510-00025	25	18	40 x 96	24
2510-00050	50	18	50 x 105	24
2510-00125	125	25	66 x 136	12
2510-00250	250	38	80 x 160	12
2510-00300	300	38	88 x 170	12
2510-00500	500	38	100 x 205	12
2510-01000	1000	38	128 x 235	6
2510-02000	2000	38	158 x 295	3
2510-03000	3000	38	184 x 340	3
2510-04000	4000	38	205 x 370	3
2510-06000	6000	38	235 x 410	3




DeLong® Style Neck without Stainless Steel Closures, without Baffles

Cat. No.	Capacity (ml)	Closures (mm)	Ø x Height (mm)	Cs.
2511-00025	25	18	40 x 96	24
2511-00050	50	18	50 x 105	24
2511-00125	125	25	66 x 136	12
2511-00250	250	38	80 x 160	12
2511-00300	300	38	88 x 170	12
2511-00500	500	38	100 x 205	12
2511-01000	1000	38	128 x 235	6
2511-02000	2000	38	158 x 295	3
2511-04000	4000	38	205 x 370	3
2511-06000	6000	38	235 x 410	3



**A large variety of caps and closures
available on request
(KAP-UTS™, stainless steel, silicone/sponge)**



Price: Prices in Euro per unit, excluding VAT.

Freight costs: Germany: We charge Euro 20.00 for delivery. Large delivery on a palette: please enquire.
Below a net order value of Euro 100.00 an additional charge of Euro 8.00 will apply.
Other countries: Please contact us for further information.

Payment: 30 days net or against pro forma invoice.

Prices are subject to change without notice. Prices supersede all prices in previous documents.
No responsibility is taken for the accuracy of the information.

01/01/2018