

Erlenmeyer-, Fernbach- and Shake Flasks from Borosilicate Glass

- complete price list on request -

Erlenmeyer Flasks (graduated)

For general laboratory use such as media preparation, storage and cultivation of microorganisms

Narrow Mouth with Break-Resistant Lip

Cat. No.	Capacity (ml)	Stopper Size	Ø x Height (mm)	Cs.
2504-00025	25	0	39 x 66	24
2504-00050	50	1	50 x 80	24
2504-00125	125	5	66 x 115	12
2504-00250	250	6	80 x 140	12
2504-00300	300	6	88 x 148	12
2504-00500	500	7	100 x 183	12
2504-01000	1000	9	128 x 223	6
2504-02000	2000	10	158 x 290	3
2504-03000	3000	10	184 x 305	3
2504-04000	4000	10	205 x 355	3
2504-06000	6000	10	237 x 395	3



Narrow Mouth

Cat. No.	Capacity (ml)	Stopper Size	Ø x Height (mm)	Cs.
2505-00125	125	5	66 x 107	12
2505-00250	250	6	80 x 133	12
2505-00500	500	7	100 x 173	6
2505-01000	1000	9	128 x 214	6



Narrow Mouth with Phenolic Rubber Lined Screw Cap

- Long-skirted screw caps protect against contamination

Cat. No.	Capacity (ml)	Cap (mm)	Ø x Height (mm)	Cs.
2514-00025	25	18-415	40 x 65	24
2514-00050	50	24-415	50 x 85	24
2514-00125	125	24-415	66 x 125	12
2514-00250	250	38-415	80 x 143	12
2514-00300	300	38-415	88 x 150	12
2514-00500	500	38-415	100 x 190	12
2514-01000	1000	38-415	128 x 220	6
2514-02000	2000	38-415	158 x 275	3
2514-03000	3000	38-415	184 x 320	3
2514-04000	4000	38-415	205 x 350	3
2514-06000	6000	38-415	235 x 390	3



Wide Mouth with GL 45 Screw Cap and Pouring Ring

- With 45 mm Polypropylene Linerless Screw Cap and Pouring Ring

Cat. No.	Capacity (ml)	Cap (mm)	Ø x Height (mm)	Cs.
2515-00250	250	45	80 x 117	12
2515-00500	500	45	100 x 156	12
2515-01000	1000	45	128 x 200	6
2515-02000	2000	45	158 x 260	3
2515-04000	4000	45	205 x 325	1
2515-06000	6000	45	235 x 370	1



Fernbach Flasks

For culturing organism requiring a large surface to volume ratio

DeLong® Style Neck, Ø 38 mm

- Without Baffles

Cat. No.	Capacity (ml)	Ø x Height (mm)	Cs.
2550-01800	1800	200 x 175	3
2550-02800	2800	205 x 250	3



Wide Mouth with Break-Resistant Lip

- With three Side Baffles

Cat. No.	Capacity (ml) and Description	Ø x Height (mm)	Cs.
2551-01800	1800, 3 Side Baffles	200 x 157	3
2551-02800	2800, 3 Side Baffles	210 x 225	3



DeLong® Style Neck, Ø 38 mm

- With three Bottom Baffles, or three Bottom and three Side Baffles

Cat. No.	Capacity (ml) and Description	Ø x Height (mm)	Cs.
2554-01800	1800, 3 Bottom Baffles	200 x 175	3
2554-02800	2800, 3 Bottom Baffles	210 x 250	3
2554-02833	2800, 6 Baffles	200 x 250	3



Wide Mouth with GL 45 Screw Cap and Pouring Ring

- 45 mm Polypropylene Screw Cap
- With three Bottom Baffles, or three Bottom and three Side Baffles

Cat. No.	Capacity (ml) and Description	Ø x Height (mm)	Cs.
2552-01800	1800, 3 Bottom Baffles	200 x 157	1
2552-02800	2800, 3 Bottom Baffles	210 x 225	1
2552-02833	2800, 6 Baffles	210 x 225	1



Wide Mouth with GL 45 Membrane Screw Cap and Pouring Ring

- 45 mm Polypropylene Screw Cap with 0,2 µm PTFE Membrane
- With three Bottom Baffles, or three Bottom and three Side Baffles

Cat. No.	Capacity (ml) and Description	Ø x Height (mm)	Cs.
2552-18000	1800, no Baffles	200 x 157	1
2552-28000	2800, no Baffles	210 x 225	1
2552-01803	1800, 3 Bottom Baffles	200 x 157	1
2552-02803	2800, 3 Bottom Baffles	210 x 225	1
2552-62180	1800 ml, 6 Baffles	200 x 157	1
2552-62280	1800 ml, 6 Baffles	210 x 225	1
5637-00046	GL 45, Polypropylene, 0,2 µm PTFE Membrane Cap		5

Shake Flasks

For preliminary fermentation and bacterial growth studies

DeLong® Style Neck, three Side Baffles

Cat. No.	Capacity (ml)	Closures (mm)	Ø x Height (mm)	Cs.
2540-00050	50	18	50 x 105	12
2540-00125	125	25	66 x 136	12
2540-00250	250	38	80 x 160	12
2540-00300	300	38	88 x 170	12
2540-00500	500	38	100 x 205	12
2540-01000	1000	38	128 x 235	6
2540-02000	2000	38	158 x 290	3



Break-Resistant Lip, three Side Baffles

Cat. No.	Capacity (ml)	Ø x Height (mm)	Cs.
2541-00125	125	66 x 116	12
2541-00250	250	80 x 140	12
2541-01000	1000	128 x 220	6
2541-02000	2000	158 x 290	3



Break-Resistant Lip, three Extra Deep Side Baffles

Cat. No.	Capacity (ml)	Ø x Height (mm)	Cs.
2542-00250	250	80 x 140	12
2542-00300	300	88 x 150	12
2542-00500	500	100 x 183	12
2542-01000	1000	128 x 220	6
2542-02000	2000	158 x 290	3



DeLong® Style Neck, three Extra Deep Side Baffles


Cat. No.	Capacity (ml)	Closures (mm)	Ø x Height (mm)	Cs.
2543-00050	50	18	50 x 105	12
2543-00125	125	25	66 x 136	12
2543-00250	250	38	80 x 160	12
2543-00300	300	38	88 x 170	12
2543-00500	500	38	100 x 205	12
2543-01000	1000	38	128 x 235	6
2543-02000	2000	38	158 x 295	3
2543-03000	3000	38	184 x 340	3
2543-04000	4000	38	205 x 370	3



DeLong® Style Neck, six Baffles

- With three Bottom Baffles, or three Bottom and three Side Baffles

Cat. No.	Capacity (ml)	Closures (mm)	Ø x Height (mm)	Cs.
2543-60250	250	38	80 x 160	12
2543-60500	500	38	100 x 205	12
2543-61000	1000	38	128 x 235	6
2543-62000	2000	38	158 x 295	3
2543-63000	3000	38	184 x 340	3
2543-64000	4000	38	205 x 370	3



Wide Mouth with GL 45 PP Membrane Screw Cap and Pouring Ring

- 45 mm Polypropylene Screw Cap with 0,2 µm PTFE Membrane (5637-00046)
- With three Bottom Baffles, or three Bottom and three Side Baffles

Cat. No.	Capacity (ml)	Closures (mm)	Ø x Height (mm)	Cs.
2545-00250	250	45	80 x 117	12
2545-00500	500	45	100 x 156	12
2545-01000	1000	45	128 x 200	6
2545-02000	2000	45	158 x 260	3



DeLong® Style Neck, four Extra Deep Bottom Baffles, with Constriction

Cat. No.	Capacity (ml)	Closures (mm)	Ø x Height (mm)	Cs.
2547-00250	250	38	81 x 160	12
2547-00500	500	38	100 x 203	12
2547-01000	1000	38	128 x 235	6



DeLong® Style Neck with Stainless Steel Closures, without Baffles

Cat. No.	Capacity (ml)	Closures (mm)	Ø x Height (mm)	Cs.
2510-00025	25	18	40 x 96	24
2510-00050	50	18	50 x 105	24
2510-00125	125	25	66 x 136	12
2510-00250	250	38	80 x 160	12
2510-00300	300	38	88 x 170	12
2510-00500	500	38	100 x 205	12
2510-01000	1000	38	128 x 235	6
2510-02000	2000	38	158 x 295	3
2510-03000	3000	38	184 x 340	3
2510-04000	4000	38	205 x 370	3
2510-06000	6000	38	235 x 410	3




DeLong® Style Neck without Stainless Steel Closures, without Baffles

Cat. No.	Capacity (ml)	Closures (mm)	Ø x Height (mm)	Cs.
2511-00025	25	18	40 x 96	24
2511-00050	50	18	50 x 105	24
2511-00125	125	25	66 x 136	12
2511-00250	250	38	80 x 160	12
2511-00300	300	38	88 x 170	12
2511-00500	500	38	100 x 205	12
2511-01000	1000	38	128 x 235	6
2511-02000	2000	38	158 x 295	3
2511-04000	4000	38	205 x 370	3
2511-06000	6000	38	235 x 410	3



A large variety of caps and closures
available on request
(KAP-UTS™, stainless steel, silicone/sponge)



01.01.2019