

# Caps and Closures

- complete price list on request -

## KAP-UTS™ Disposable Closures

- For culture tubes and flasks with DeLong® style necks and access ports
- From nontoxic autoclavable Polypropylene
- Seven colours for easy identification
- Six sizes from 13 mm to 38 mm available
- Same evaporation and sterility characteristics as stainless steel closures
- Molded pressure fins for firm grip
- Autoclavable to 140 °C




Cat. No.	Size	Colour	Cs.
2007-13001	13 mm	black	1000
2007-13002	13 mm	red	1000
2007-13003	13 mm	yellow	1000
2007-13004	13 mm	blue	1000
2007-13005	13 mm	white	1000
2007-13006	13 mm	green	1000
2007-16001	16 mm	black	1000
2007-16002	16 mm	red	1000
2007-16003	16 mm	yellow	1000
2007-16004	16 mm	blue	1000
2007-16005	16 mm	white	1000
2007-16006	16 mm	green	1000
2007-18001	18 mm	black	1000
2007-18002	18 mm	red	1000
2007-18003	18 mm	yellow	1000
2007-18004	18 mm	blue	1000
2007-18005	18 mm	white	1000
2007-18006	18 mm	green	1000
2007-20001	20 mm	black	500
2007-20002	20 mm	red	500
2007-20003	20 mm	yellow	500
2007-20004	20 mm	blue	500
2007-20005	20 mm	white	500
2007-20006	20 mm	green	500
2007-25001	25 mm	black	500
2007-25002	25 mm	red	500
2007-25003	25 mm	yellow	500
2007-25004	25 mm	blue	500
2007-25005	25 mm	white	500
2007-25006	25 mm	green	500
2007-25013	25 mm	colourless	500
2007-38005	38 mm	white	250

## Stainless Steel Closures

- For aerobic culture in tubes and flasks utilizing DeLong® Necks
- Made from 305 stainless steel (prevents corrosion)
- With finger grips for firm fit on the flask
- Suitable for dry heat sterilization and repeated autoclaving

Cat. No.	Tube outside-Ø	Cs.
2005-00013	13 mm	144
2005-00016	16 mm	144
2005-00018	18 mm	144
2005-01824	18 mm	24
2005-00020	20 mm	144
2005-00025	25 mm	72
2005-02512	25 mm	12
2005-00038	38 mm	36
2005-03803	38 mm	3
2005-03806	38 mm	6
2005-03812	38 mm	12



## Silicone/Sponge Closures

- Sealed-in silicone sponge filter is ideal for aerobic cultures
- Flanged to fit snugly over outside of vessel neck preventing passage of aerosols
- Closures will withstand autoclaving and dry heat sterilization (180 °C)

Cat. No.	Description	Size (mm)	Cs.
2004-00002	13 - 18 mm Culture tubes 18 mm screw necks	13-18	10
2004-00003	18 - 25 mm Culture tubes 25 - 50 ml Erlenmeyer flasks 25 - 125 ml DeLong® flasks	18-25	10
2004-00004	125 - 300 ml Erlenmeyer flasks	28	10
2004-00005	38 mm Culture tubes 500 - 1000 ml Erlenmeyer flasks 250 - 6000 ml DeLong® flasks Fit over 45 mm Screw Necks	38	10



## Silicone/Sponge Plug Closures

- Tapered silicone sponge with smooth exterior for tight seal
- Extra long for easy insertion and removal

Cat. No.	Description	Size (mm)	Cs.
1924-00032	For 25 mm culture tubes	32	25



Find more products for your lab from:



BELLCO



Sterilin Limited



SciLabware



porvair sciences



EXCEL Scientific



Life Sciences

MBL® PYREX® QUICKFIT® AZLON® Polystyrol (PS), Polypropylen (PP), PTFE, HDPE, LDPE, PFA, PVC, PC,



