

Laboratory equipment suitable for gel and blot applications

Shaker "Belly Dancer"



- Suitable for gel and blot applications
- Continuous gentle or thorough shaking
- Unique three dimensional motion guarantees accurate mixing action and offers the possibility to stack blotting trays
- Durable construction, low noise level
- Ambient temperature operation: 0 °C to 40 °C
- Can be used in CO₂ incubators and cool rooms
- Speed control: 0 - 100 rpm
- Maximum load: 9.1 kg
- Platform size: 300 x 300 mm
- Dimensions (w x d x h): 381 x 381 x 229 mm
- Weight: 7.7 kg
- Additional platform can be added (see picture)

Cat.-No.	Description
BDRLS0006	"Belly Dancer", with finger mat for tubes
BDRLS0004	"Belly Dancer", with additional platform (approximately 300 x 300 mm) (see picture)
BDRLS0005	"Belly Dancer", with extended platform (approximately 450 x 450 mm)

Shaker "Belly Button"



- Suitable for gel and blot applications
- Very compact version of the "Belly Dancer" (see above), fits also into smaller incubators
- Platform tilt is continuously adjustable from level to 8 °
- Ambient temperature operation: 0 °C to 60 °C
- Can be used in CO₂ incubators and cool rooms
- Speed control: 6 - 100 rpm
- Maximum load: 2.3 kg
- Dimensions (w x d x h): 203 x 203 x 133 mm
- Optionally available: non-skid mat for shaking tubes (see picture)

Cat.-No.	Description
BBUAAUV1S	"Belly Button"
ACCBB0001	Non-skid mat for tubes, 178 x 178 mm, black

Roller Mixers



- Space saving low profile design
- Available with 6 or 11 rolls (6 rolls model can be extended to 11 rolls)
- Suitable for culture and hybridization bottles or tubes due to removable rolls and CO₂ incubator, dimension of one roll (Ø x L): 29 x 254 mm
- Can be used with 3 big culture bottles (850 cm²), 4 hybridization bottles or with up to 30 tubes (10 x 75 mm), if used with 11 rolls
- Speed control: 2 - 38 rpm
- Maximum load: 4 liter liquid
- Ambient temperature operation: 0 °C to 60 °C, can be used in CO₂ incubator
- Dimensions (w x d x h): 362 x 381 x 60.5 mm

Cat.-No.	Description
ROLAAUV1S	Roller mixer with 6 rolls
ROLLS0001	Roller mixer with 11 rolls
ACCRLSR10	Single roll


Shaking Water Bath



- Combination of modified "Belly Dancer" (see front page) and hybridization water bath
- Units can also be used separately
- Separation of shaker and water bath offers the possibility for easy cleaning and decontamination when working with radioactive material
- The PET clear lid has an easy grip handle and a thermometer port
- PID microprocessor heat control, safety thermostat
- Temperature range: 30 °C to 75 °C
- Quick heating, 3 °C rise per minute
- Excellent temperature control and stability: ± 0.2 °C

Cat.-No.	Description
BWBAA220S	Shaking Water Bath

Controller "Blot Washer"




Belly Dancer
(has to be ordered separately)

- The "Blot Washer" is a controller with integrated pump, ideal for the automation of multiple washing steps for Western Blots
- Up to 20 washing programs can be created and saved
- Programmable washing time from seconds to 99 hours
- Simultaneous processing of up to four containers with different sizes
- Liquid dispensing rate: approx. 2 ml/second
- Liquid aspirating rate: approx. 5 ml/second
- Delivered with 2 x 4 liter bottles, 2 blot trays and vacuum tube with connectors
- Can be combined with the "Belly Dancer" (see front page, not included, has to be ordered separately)
- Accessories for "Blot-Washer" (tubing, clamps, clips, trays, etc.) on request

Cat.-No.	Description
BLWAA220S	"Blot Washer" , controller with integrated pump

"Hi/Lo Profile Rocker"



- Space saving low profile design, stainless steel
- Autoclavable (without motor)
- Gentle to vigorous shaking motion
- Speed control: 0 - 22 rpm
- Platform tilt is continuously adjustable from: 2 ° - 11 °
- Ambient temperature operation: 0 °C to 60 °C
- Can be used in CO₂ incubators and cool rooms,
- Maximum load: 9.1 kg, max. operating humidity: 95%
- Platform with non-skid mat: 305 x 305 mm
- Option: extended capacity with up to 3 additional platforms
- Distance between platforms (if more than one is used): 114 mm
- Dimensions (w x d x h): 349 x 400 x 44.5 mm (without additional platforms)

Cat.-No.	Description
ROCAA220S	"Hi/Lo Profile Rocker"
ROCLS0004	"Hi/Lo Profile Rocker" , with 1 additional platform
ROCLS0005	"Hi/Lo Profile Rocker" , with 2 additional platforms
ROCLS0006	"Hi/Lo Profile Rocker" , with 3 additional platforms

Please note that in Europe all equipment listed above is not delivered with a flat cable!

Price: Prices in Euro per unit, excluding VAT.
Freight costs: Germany: We charge Euro 20.00 for delivery. Large delivery on a palette: please enquire.
 Below a net order value of Euro 100.00 an additional charge of Euro 8.00 will apply.
 Other countries: Please contact us for further information.

Payment: 30 days net or against pro forma invoice.

Prices are subject to change without notice. Prices supersede all prices in previous documents. No responsibility is taken for the accuracy of the information.