

## Pipette Plugger and Serological Glass Pipettes

Plug your glass pipettes with the pipette plugger from Bellco Glass:  
Plugged serological glass pipettes for your cell culture application are an economic and ecologically sensitive solution vs. disposable plastic pipettes.  
The machine can plug 40, 50 or even 60 pipettes per minute. It has different operating modes for continuous or on-demand supply of cotton. The uniform plug shape guarantees continuously and consistently plugged pipettes. The cotton supplied with the machine will last for approximately 35,000 pipettes. The compact design of the machine and a handle make it easy to transport within your lab.



- Three plugging speeds:  
40, 50 or 60 pipettes per minute
- Operates in continuous or on-demand mode
- Uniform plug shape insures proper plug fit
- Supplied with 1 bag of cotton roving  
(Cat. No. 7750-22222 see below)

Cat. No.	Description
7750-11220	Pipette plugger (W x D x H): 25 x 28 x 24 cm
A515-A18	Scissor Blade Assembly only

### Cotton Roving for Pipette Plugger



- Absorbent, 12 grain premium long-fiber cotton
- Packed coiled in plastic bag, approx. 0,9 kg
- Plugs approximately 35.000 pipettes

Cat. No.	Description	Unit
7750-22222	Cotton Roving	approx. 0,9 kg

# Glass Pipettes (excerpt)

## Serological Pipettes with Cotton Plugging Top

- Type 1 Class B Borosilicate Glass
- Constricted top for cotton plugging (supplied without cotton)
- Tempered tip, to deliver, blow out
- Calibrated to tip

### Standard



Cat. No.	Size (ml)	Tolerance $\pm$ (ml)	Length (mm)	Colour Code	Cs.
1210-00201	0,2 x 0,01	0,008	330	Black	18
1210-01010	1,0 x 0,10	0,020	350	Red	18
1210-01001	1,0 x 0,01	0,020	350	Yellow	18
1210-02010	2,0 x 0,10	0,020	350	Green	18
1210-02001	2,0 x 0,01	0,020	350	White	18
1210-05010	5,0 x 0,10	0,040	350	Blue	18
1210-10010	10 x 0,10	0,060	350	Orange	18
1210-25010	25 x 0,10	0,100	400	White	18
1210-50020	50 x 0,20	0,400	685	Red	12

### Large Orifice



- 3 mm orifice for handling tissue homogenates and cell suspensions

Cat. No.	Size (ml)	Tolerance $\pm$ (ml)	Length (mm)	Colour Code	Cs.
1232-25010	25 x 0,1	0,1	400	White	18
1232-50020	50 x 0,2	0,5	470	Red	12

### “Shorties”

- Shorter length facilitates use in hoods



Cat. No.	Size (ml)	Tolerance $\pm$ (ml)	Length (mm)	Colour Code	Cs.
1226-01010	1 x 0,1	0,020	175	Red	18
1226-02010	2 x 0,1	0,025	175	Green	18
1226-05010	5 x 0,1	0,080	175	Blue	18
1226-25050	25 x 0,5	0,250	230	White	18

### “Shorties” with Large Orifice

- 3 mm orifice for handling tissue homogenates and cell suspensions
- Shorter length facilitates use in hoods



Cat. No.	Size (ml)	Tolerance $\pm$ (ml)	Length (mm)	Colour Code	Cs.
1228-05010	5 x 0,1	0,14	175	Blue	18
1228-25050	25 x 0,5	0,25	225	White	18

**Information on further pipettes and pipette containers upon request**