


Microscope Slide Incubation Trays

for Immunology, Histology and Molecular Biology

Price List 2019

The incubation trays can be used as staining or humidity chamber in Immunohistochemistry, Histology, *in situ* Hybridization, Microbiology, Haematology and Molecular Biology. Due to anti-slip strips, the slides are positioned firmly in the trays. The boxes can be used in refrigerators and incubators (temperature range: -20 °C to +70 °C).

Available in several sizes for different applications



**Incubation Tray
SlideShow30**

- For up to 30 slides,
slide dimensions: 26 mm x 76 mm
- Tray dimensions (w x d x h):
380 mm x 330 mm x 40 mm



**Incubation Tray
SlideShow20**

- For up to 20 slides,
slide dimensions: 26 mm x 76mm
- Tray dimensions (w x d x h):
380 mm x 230 mm x 40 mm



**Incubation Tray
SlideShow10**

- For up to 10 slides,
slide dimensions: 26 mm x 76 mm
- Tray dimensions (w x d x h):
380 mm x 150 mm x 40 mm


Features:

- Different sizes
- Special polymer rails grip slides securely
- Easy handling of the slides
- **Clear or black lids** for fluorescence work
- Wells between the rails can be filled with liquid to maintain humidity in the chamber
- Neat, stackable design to maximize storage space
- Proven reliability and long life
- Slide rails spaced 35 mm apart
- Material:
 - Tray: black ABS (Acrylonitrile Butadiene Styrene), 5 mm
 - Black lid: ABS (Acrylonitrile Butadiene Styrene), 5 mm
 - Clear lid: PETG (Polyethylene Terephthalate Glycol), 5mm
- Operation temperature range: -20 °C to +70 °C


Do not use Chlorinated Hydrocarbons or aggressive solvents!

Schematic representation of the different trays


Tray for 10 slides



Tray for 20 slides



Tray for 30 slides



| Cat. No. | Description |
|--------------|---|
| 074/4300MS10 | Incubation tray with clear lid for up to 10 slides |
| 074/4300MS11 | Incubation tray with black lid for up to 10 slides |
| 074/4300MS20 | Incubation tray with clear lid, for up to 20 slides |
| 074/4300MS21 | Incubation tray with black lid, for up to 20 slides |
| 074/4300MS30 | Incubation tray with clear lid, for up to 30 slides |
| 074/4300MS31 | Incubation tray with black lid, for up to 30 slides |

→ Replacement lids available on request!



**Dunn Labortechnik GmbH also offers
Multispot Microscope Slides for immunofluorescence,
cell culture or *in situ* hybridization**



and

**a wide range of antibodies, antigens, proteins, ELISA kits, plasma,
sera and substrates from various manufacturers:**



- Price:** Prices in Euro per piece, excluding VAT.
- Freight charges:** Germany: We charge Euro 20.00 for delivery. Large delivery on a palette please enquire. Below a net order value of Euro 100.00 an additional charge of Euro 8.00 will apply. Other countries: please contact us for further information.
- Payment terms:** 30 days net or against pro forma invoice.

Prices are subject to change without notice. Prices supersede all prices in previous documents. No responsibility is taken for the accuracy of the information.