


Microscope Slide Incubation Trays

for Immunology, Histology and Molecular Biology

Product Information 2018

The incubation trays can be used as staining or humidity chamber in Immunohistochemistry, Histology, *in situ* Hybridization, Microbiology, Haematology and Molecular Biology. Due to anti-slip strips, the slides are positioned firmly in the trays. The boxes can be used in refrigerators and incubators (temperature range: -20 °C to +70 °C).

Available in several sizes for different applications



**Incubation Tray
SlideShow30**

- For up to 30 slides,
slide dimensions: 26 mm x 76 mm
- Tray dimensions (w x d x h):
385 mm x 330 mm x 40 mm



**Incubation Tray
SlideShow20**

- For up to 20 slides,
slide dimensions: 26 mm x 76mm
- Tray dimensions (w x d x h):
385 mm x 245 mm x 40 mm



**Incubation Tray
SlideShow10**

- For up to 10 slides,
slide dimensions: 26 mm x 76 mm
- Tray dimensions (w x d x h):
385 mm x 160 mm x 40 mm


Features:

- Different sizes
- Special polymer rails grip slides securely
- Easy handling of the slides
- **Clear or black lids** for fluorescence work
- Wells between the rails can be filled with liquid to maintain humidity in the chamber
- Neat, stackable design to maximize storage space
- Proven reliability and long life
- Slide rails spaced 35 mm apart
- Material:
 - Tray: black ABS (Acrylonitrile Butadiene Styrene), 5 mm
 - Black lid: ABS (Acrylonitrile Butadiene Styrene), 5 mm
 - Clear lid: PETG (Polyethylene Terephthalate Glycol), 5mm
- Operation temperature range: -20 °C to +70 °C


Do not use Chlorinated Hydrocarbons or aggressive solvents!

Schematic representation of the different trays


Tray for 10 slides



Tray for 20 slides




Tray for 30 slides



Cat. No.	Description
074/4300MS10	Incubation tray with clear lid for up to 10 slides
074/4300MS11	Incubation tray with black lid for up to 10 slides
074/4300MS20	Incubation tray with clear lid, for up to 20 slides
074/4300MS21	Incubation tray with black lid, for up to 20 slides
074/4300MS30	Incubation tray with clear lid, for up to 30 slides
074/4300MS31	Incubation tray with black lid, for up to 30 slides

→ Replacement lids available on request!



**Dunn Labortechnik GmbH also offers
Multispot Microscope Slides for immunofluorescence,
cell culture or *in situ* hybridization**



and

**a wide range of antibodies, antigens, proteins, ELISA kits, plasma,
sera and substrates from various manufacturers:**

