


Products for Cryo Storage - Product List 2019

Expell Cryo Tubes


DNase-/RNase-free. Manufactured from FDA-approved, medical grade polypropylene.
 Manufactured in a clean room environment and according to ISO 9001, ISO 14001, ISO 13485.
 Sterilized according to ISO 11137-2:2013 – Sterility Assurance Level (SAL) by standard products 10-6.



Cat. No.	Description
5040500C	Expell cryo tubes, 0,5 ml, sterile , length: 28.6 mm, outer diameter 12.08 mm, 10 x 100 tubes
5041000C	Expell cryo tubes, 1,0 ml, sterile , length: 36.1 mm, outer diameter 12.08 mm, 10 x 100 tubes
5042000C	Expell cryo tubes, 2,0 ml, sterile , length: 47.6 mm, outer diameter 12.08 mm, 10 x 100 tubes
5042001C	Expell cryo tubes, 2,0 ml , non-sterile, length: 47.6 mm, outer diameter 12.08 mm, 2 x 500 tubes
5044000C	Expell cryo tubes, 4,0 ml, sterile , length: 73.3 mm, outer diameter 12.08 mm, 5 x 100 tubes
5045000C	Expell cryo tubes, 5,0 ml, sterile , length: 88.5 mm, outer diameter 12.08 mm, 5 x 100 tubes
504W000C	Coloured insert for Expell cryo tubes, white, 500 pcs.
504G000C	Coloured insert for Expell cryo tubes, green, 500 pcs.
504R000C	Coloured insert for Expell cryo tubes, red, 500 pcs.
504B000C	Coloured insert for Expell cryo tubes, blue, 500 pcs.
504V000C	Coloured insert for Expell cryo tubes, purple, 500 pcs.
504Y000C	Coloured insert for Expell cryo tubes, yellow, 500 pcs.

Expell Cryo Boxes

- Made of high quality durable polypropylene, suitable for storing at -90 °C and autoclaving at 121 °C for 20 minutes.
- Numerical coding for each tube for quick and easy sample identification.
- Star-foot lock at the bottom of each box well, enabling one-handed tube opening without taking the tube out of the box.
- Slanted edges for easier identification and handling.
- Save up to 23 % storage space in combination with the Expell cryo tubes. In addition, use of 0.5 ml tubes allows stacking of the boxes.



Cat. No.	Description	Case
5040510C	Expell cryo box, 10 x 10 positions, for 0.5 ml Expell cryo tubes, box dimensions (w x l x h): 130 mm x 130 mm x 32 mm	180 boxes
5041010C	Expell cryo box, 10 x 10 positions, for 1 ml Expell cryo tubes, box dimensions (w x l x h): 130 mm x 130 mm x 40 mm	105 boxes
5042010C	Expell cryo box, 10 x 10 positions, for 2 ml Expell cryo tubes, box dimensions (w x l x h): 130 mm x 130 mm x 50 mm	105 boxes
5044010C	Expell cryo box, 10 x 10 positions, for 4 ml Expell cryo tubes, box dimensions (w x l x h): 130 mm x 130 mm x 80 mm	90 boxes
5045010C	Expell cryo box, 10 x 10 positions, for 5 ml Expell cryo tubes, box dimensions (w x l x h): 130 mm x 130 mm x 90 mm	60 boxes

Capp Liquid handling products and Rondo laboratory equipment. Please contact us for further information.



Price: Prices in Euro per unit, excluding VAT.

Freight charges: Germany: We charge Euro 20.00 for delivery. Large delivery on a palette please enquire.
Below a net order value of Euro 100.00 an additional charge of Euro 8.00 will apply.
Other countries: please contact us for further information.

Payment terms: 30 days net or against pro forma invoice.

Prices are subject to change without notice. Prices supersede all prices in previous documents. No responsibility is taken for the accuracy of the information.