

“Intelli-Plates” for protein crystallography applications

- Molecular structural analysis
- Three-dimensional structure of proteins
- Structure based drug research

96-well INTELLI-PLATE®		
 <p>Catalog #: 102-0011-00 Original INTELLI-PLATE Well Volume: 4 µl and 10 µl Reservoir volume: 300 µl</p>	 <p>Catalog #: 102-0001-00 96-2 Low volume reservoir Well Volume: 4 µl and 10 µl Reservoir Volume: 100 µl</p>	 <p>Catalog #: 102-0001-10 96-2 Low profile Well Volume: 4 µl and 10 µl Reservoir Volume: 100 µl</p>
 <p>Catalog #: 102-0001-03* 96-3 Low volume reservoir Well Volume: 1 µl Reservoir Volume: 100 µl</p>	 <p>Catalog #: 102-0001-13* 96-3 low profile INTELLI-PLATE Well Volume: 1 µl Reservoir Volume: 100 µl</p>	 <p>Catalog #: 102-0001-01 Flat bottom Well Volume: N/A Reservoir Volume: 300 µl</p>
<p>UV Compatible Plates*</p> <p>102-0001-03 102-0001-13 102-0001-20 102-0001-21</p>		
 <p>Catalog #: 102-0001-20* 96-2 Shallow well Well Volume: 2 µl Reservoir Volume 100 µl</p>	 <p>Catalog #: 102-0001-21* 96-2 Shallow well low profile Well Volume: 2 µl Reservoir Volume: 100 µl</p>	
24-well INTELLI-PLATE®		48-well INTELLI-PLATE®
 <p>Catalog #: 102-0004-00 24-4 Well Volume: 5 µl Reservoir Volume: 650 µl</p>	 <p>Catalog #: 102-0003-00 48-3 Screening Well Volume: 4 µl Reservoir Volume: 500 µl</p>	 <p>Catalog #: 102-0002-00 48-2 Optimization Well Volume: 4 µl and 20 µl Reservoir Volume: 500 µl</p>

- **96-well “Intelli-Plate” (8 versions):**
 - two crystallization wells (max. 4 µl and 10 µl) and optionally 100 µl or 300 µl reservoir
 - three crystallization wells (1 µl) and 100 µl reservoir
 - low profile and shallow well plates, now UV-compatible
 - highly polished, flat crystallization area and 300 µl reservoir
- **48-well “Intelli-Plate” (2 versions):**
 - for optimization with two wells (4 µl and 20 µl) and 500 µl reservoir
 - for screening with three wells (4 µl) and 500 µl reservoir
- **24-well “Intelli-Plate”:**
 - for “sitting drop”, with four wells (5 µl) and 650 µl reservoir



- All plates available with barcode
- Standard SBS microplate footprint
- Single packaged to prevent electrostatic charging of the plates

We also offer compatible sealing foils (ThermalSeal RTS) from Excel Scientific to protect your protein samples. Please contact us or visit our webpage for further information.

“Intelli-Plates” Price Information

Cat. No.	Description	Qty.
102-0011-00	96-well Original “Intelli-Plate”, PS, Reservoir 300 µl, 2 Protein Wells max. 10 µl and 4 µl	12 x 10
102-0011-00S	dto., small pack	1 x 10
102-0011-00B	dto., bulk pack	120 x 10
102-0001-00	96-well “Intelli-Plate”, LBR, Reservoir 100 µl, 2 Protein Wells max. 10 µl and 4 µl	12 x 10
102-0001-00S	dto., small pack	1 x 10
102-0001-00B	dto., bulk pack	120 x 10
102-0001-10	96-well “Intelli-Plate”, Low Profile , LBR, Reservoir 100 µl, 2 Wells max. 10 µl and 4 µl	6 x 20
102-0001-10S	dto., small pack	1 x 20
102-0001-10B	dto., bulk pack	60 x 20
102-0001-03	96-well “Intelli-Plate”, LBR, Reservoir 100 µl, 3 Wells à 1 µl, <u>UV compatible</u>	12 x 10
102-0001-03S	dto., small pack	1 x 10
102-0001-03B	dto., bulk pack	120 x 10
102-0001-13	96-well “Intelli-Plate”, Low Profile , LBR, Reservoir 100 µl, 3 Wells à 1 µl, <u>UV compatible</u>	6 x 20
102-0001-13S	dto., small pack	1 x 20
102-0001-13B	dto., bulk pack	60 x 20
102-0001-01	96-well “Intelli-Plate”, PS, Reservoir 300 µl, with shelf	12 x 10
102-0001-01S	dto., small pack	1 x 10
102-0001-01B	dto., bulk pack	120 x 10
102-0004-00	24-well “Intelli-Plate”, LBR, Reservoir 650 µl, 4 Wells à 5 µl	12 x 10
102-0004-00S	dto., small pack	1 x 10
102-0004-00B	dto., bulk pack	120 x 10
102-0002-00	48-well “Intelli-Plate” for optimization, PS, Reservoir 500 µl, 2 Wells à 4 µl and 20 µl	12 x 10
102-0002-00S	dto., small pack	1 x 10
102-0002-00B	dto., bulk pack	120 x 10
102-0003-00	48-well “Intelli-Plate” for screening, PS, Reservoir 500 µl, 3 Wells à 4 µl	12 x 10
102-0003-00S	dto., small pack	1 x 10
102-0003-00B	dto., bulk pack	120 x 10
102-0001-20	96-Well “Intelli-Plate”, 2 Shallow Wells à 2 µl, Reservoir 100 µl, LBR, <u>UV compatible</u>	12 x 10
102-0001-20S	dto., small pack	1 x 10
102-0001-20B	dto., bulk pack	120 x 10
102-0001-21	96-Well “Intelli-Plate”, Low Profile , 2 Shallow Wells à 2 µl, Reservoir 100 µl, LBR, <u>UV compatible</u>	6 x 20
102-0001-21S	dto., small pack	1 x 20
102-0001-21B	dto., bulk pack	60 x 20
****_****_**B-BC	All versions of the 96-well “Intelli-Plate” with bar code (except low profile plates)	240 plates
****_****_**B-BC-LP	All low profile versions of the 96-well “Intelli-Plate” with bar code	500 plates