


## Wash Bottles in many different versions





Cat. No.	Description	Cs.
	<b>Wash bottles, classic style, narrow neck, LDPE (**),</b> drip-free <b>PP cap</b> , square shoulder, possible to use precision jet tip in every position	
	<b>Volume</b> <b>Cap</b>	
WGF230P	250 ml                      white	5
WGF236P	500 ml                      white	5



Cat. No.	Description	Cs.
	<b>Wash bottles, narrow neck, sloped shoulder, LDPE (**),</b> drip-free, possible to use precision jet tip in every position, drip-free white <b>PP cap</b> , 3 sizes	
	<b>Volume</b> <b>Cap</b>	
WGF330P	250 ml                      white	5
WGF330VRAC	250 ml                      white	bulk package 150
WGF336P	500 ml                      white	5
WGF336VRAC	500 ml                      white	bulk package 100
WGF240P	1000 ml                      white	5
WGF240VRAC	1000 ml                      white	bulk package 55



Cat. No.	Description	Cs.
	<b>Wash bottles, printed with a choice of 7 popular</b> <b>laboratory solvents, and appropriate hazard symbols,</b> <b>LDPE (**),</b> possible to use precision jet tip in every position, available in 2 sizes, 6 colour coded drip-free <b>PP caps</b>	
	<b>Volume</b> <b>Cap</b> <b>Label</b>	
WGF500	250 ml                      white	-                      5
WGF501	500 ml                      white	-                      5
WGF530P	250 ml                      white	"Water"                      5
WGF531P	250 ml                      red	"Acetone"                      5
WGF532P	250 ml                      green	"Methanol"                      5
WGF533P	250 ml                      blue	"Isopropanol"                      5
WGF534P	250 ml                      yellow	"I.M.S."                      5
WGF535P	250 ml                      orange	"Ethanol"                      5
WGF536P	500 ml                      white	"Water"                      5
WGF537P	500 ml                      red	"Acetone"                      5
WGF538P	500 ml                      green	"Methanol"                      5
WGF539P	500 ml                      blue	"Isopropanol"                      5
WGF540P	500 ml                      yellow	"I.M.S."                      5
WGF541P	500 ml                      orange	"Ethanol"                      5
WGF600P	250 ml                      white	"Dist. Water"                      5
WGF601P	500 ml                      white	"Dist. Water"                      5



Cat. No.	Description	Cs.
	<b>Wash bottles, coloured, narrow neck, LDPE (**),</b> <b>drip-free PP caps</b> , possible to use precision jet tip in every position, available in 2 colours, tip has the same colour as the bottle	
	<b>Volume</b> <b>Bottle</b>	
WGR038P	500 ml                      red	5
WGR138P	500 ml                      blue	5



Cat. No.	Description	Cs.
	<b>Wash bottles, non-venting, wide neck, LDPE (**),</b> with blank writing space, printed in 3 languages, supplied with adhesive hazard labels	
WGW001PML	250 ml	5
WGW002PML	500 ml	5
WGW020PML	1000 ml	5
WGWPEN01	Marker pen, black, for wash bottle	2



Cat. No.	Description	Cs.
	<b>Wash bottles, patented DripLok® vapour venting technology, wide neck, LDPE (**), precision jet tip aperture with white PP cap, with blank writing space</b>	
WGW001VTML	250 ml	5
WGW002VTML	500 ml	5
WGW004VTML	500 ml, Caps in different colours	5
WGW020VTML	1000 ml	5
WGWPEN01	Marker pen, black, for wash bottle	2



	<b>Wash bottles, patented DripLok® vapour venting technology wide neck, LDPE (**), labelled for gram staining</b>	
WGW0255V	250 ml Bottles with label: "Crystal violet", "Lugol's iodine", "Gram iodine", "Gram colour/ de-colourising" and "Safranin solution"	Set/5



	<b>Wash bottles, wide neck, LPDE (**), colour coded leak proved cap (5 colours), precision jet tip aperture and PP Cap, 3 bottle sizes</b>		
	<b>Volume</b>	<b>Cap</b>	
WGW230P	250 ml	white	5
WGW231P	250 ml	red	5
WGW232P	250 ml	yellow	5
WGW233P	250 ml	green	5
WGW234P	250 ml	blue	5
WGW235P	500 ml	blue	5
WGW236P	500 ml	white	5
WGW237P	500 ml	red	5
WGW238P	500 ml	yellow	5
WGW239P	500 ml	green	5
WGW443P	1000 ml	white	5
WGW444P	1000 ml	yellow	5
WGW445P	1000 ml	red	5
WGW446P	1000 ml	blue	5
WGW447P	1000 ml	green	5
WGW229P	250 ml bottles, one of each colour		Set/5



	<b>Same features as above, patented DripLok® vapour venting technology</b>		
	<b>Volume</b>	<b>Cap</b>	
WGW234PVENT	500 ml	orange	5
WGW235PVENT	500 ml	blue	5
WGW236PVENT	500 ml	white	5
WGW237PVENT	500 ml	red	5
WGW238PVENT	500 ml	yellow	5
WGW239PVENT	500 ml	green	5
WGW443VT	1000 ml	white	5
WGW444VT	1000 ml	yellow	5
WGW445VT	1000 ml	red	5
WGW446VT	1000 ml	blue	5
WGW447VT	1000 ml	green	5
WGW448VT	1000 ml	orange	5

**Please ask for our complete price list  
on the Azlon product range from SciLabware!**



Cat. No.	Description	Cs.
<b>Wash bottles, LDPE (**), Label, writing space according to CPA (Clinical Pathology Accreditations), removable tip, venting cap (volatile solutions), wide neck</b>		
WTR250VT	250 ml Venting cap, red	5
WTR500VT	500 ml Venting cap, red	5



**Wash bottles, wide neck, LDPE (\*\*), label with 4 languages (English, French, German and Spanish°), non-venting, colour coded cap, precision jet tip aperture with one-piece PP-closure**

	Volume	Cap	Label	Cs.
WGW530PML	250 ml	white	"Water"	5
WGW531PML	250 ml	red	"Acetone"	5
WGW532PML	250 ml	green	"Methanol"	5
WGW533PML	250 ml	blue	"Isopropanol"	5
WGW534PML	250 ml	yellow	"I.M.S."	5
WGW535PML	250 ml	orange	"Ethanol"	5
WGW401PML	500 ml	white	"Hypochlorid"	5
WGW536PML	500 ml	white	"Water"	5
WGW537PML	500 ml	red	"Aceton"	5
WGW538PML	500 ml	green	"Methanol"	5
WGW539PML	500 ml	blue	"Isopropanol"	5
WGW540PML	500 ml	yellow	"I.M.S."	5
WGW541PML	500 ml	orange	"Ethanol"	5
WGW600PML	250 ml	white	"Distilled Water"	5
WGW601PML	500 ml	white	"Distilled Water"	5
WGW602PML	250 ml	white	"Deionized Water"	5
WGW603PML	500 ml	white	"Deionized Water"	5
WGW700PP	500 ml	white	"Azlon"	5
WGW701PP	500 ml	white	"Dichloromethane"	5
WGW703PP	500 ml	white	"Methyl ethyl ketone"	5

°other languages such as Arabic, Russian or Polish are also available



**Same as above., patented DripLok® vapour venting technology**

	Volume	Cap	Label	Cs.
WGW531VTML	250 ml	red	"Acetone"	5
WGW532VTML	250 ml	green	"Methanol"	5
WGW533VTML	250 ml	blue	"Isopropanol"	5
WGW534VTML	250 ml	yellow	"I.M.S."	5
WGW535VTML	250 ml	orange	"Ethanol"	5
WGW537VTML	500 ml	red	"Acetone"	5
WGW538VTML	500 ml	green	"Methanol"	5
WGW539VTML	500 ml	blue	"Isopropanol"	5
WGW540VTML	500 ml	yellow	"I.M.S."	5
WGW541VTML	500 ml	orange	"Ethanol"	5



**Round integral wash bottles with an integrated spout prevent contamination**

	Volume	Cap	Label	Cs.
WGW800P	500 ml	mixed	No Label	5
WGW801P	500 ml	white	"Distilled Water"	5
WGW802P	500 ml	red	"Acetone"	5
WGW803P	500 ml	green	"Methanol"	5
WGW804P	500 ml	white	"Isopropanol"	5
WGW805P	500 ml	orange	"Ethanol"	5