

## FLEXIBLE SPEED. HIGH-CAPACITY PERFORMANCE.

Supports a broad range of applications, from blood processing to viral cultures, with customizable speed settings and drum capacities to match your lab's workflow.

### MODELS

- Standard Brush-less, Low Speed
- Standard Brush-less
- Low Profile Brush-less, Low Speed
- Low Profile Brush-less

### EASY MOUNTING BRACKET

### DRUM TILT

7°

### AVAILABLE VOLTAGE

115V  
230V

### DIGITAL DISPLAY

### SPEED RANGE

2 - 60 rpm  
0.5 - 15 rpm  
Custom Speeds  
Available





Part # 7736-20025

**156 POSITION  
25mm TUBES**



Part # 7736-20050

**128 POSITION  
30mm TUBES**



Part # 7736-20164

**164 POSITION  
16mm TUBES**



Part # 7736-20351

**351 POSITION  
16mm TUBES**



Part # 7736-21106

**106 POSITION  
LOW PROFILE  
18mm TUBES**



**CONTACT US:**

**Tel.:** +49 (0) 26 83 / 4 30 94

**E-mail:** [info@dunnlab.de](mailto:info@dunnlab.de)

**Internet:** [www.dunnlab.de](http://www.dunnlab.de)

