

Special Global N-Biotek Sale

25 % Discount on All Models!

Shakers, Incubators, Shaking Incubators and More!

from August to November 2026

Incubators



Features

- Temperature range: 5 °C above ambient to 60 °C (± 0.25 °C at 37 °C)
- CO₂ range: 0 to 20 % (infrared sensor); CO₂ accuracy: ± 0.1 % at 5 %
- Stainless steel interior and perforated shelves
- Supplied with water tray for natural humidity
- Many optional accessories available such as additional shelves, access port, stacking kits, UV or hot air sterilization, Peltier cooling, oxygen control, and more!

Standard Models

Model	Chamber Volume	External Dimensions (w x d x h)	Internal Dimensions (w x d x h)	Shelves
NB203M	15.2 l	292 x 333 x 433 mm	224 x 200 x 340 mm	2 (max. 3)
NB203	42 l	408 x 482 x 550 mm	320 x 350 x 370 mm	2 (max. 4)
NB203XL	180 l	560 x 665 x 945 mm	473 x 528 x 710 mm	3 (max. 6)
NB203XXL	850 l	820 x 950 x 1840 mm	698 x 799 x 1528 mm	3 (max. 15)

Models with Touch Screen

Model	Chamber Volume	External Dimensions (w x d x h)	Internal Dimensions (w x d x h)	Shelves
NBT203	42 l	408 x 482 x 550 mm	320 x 350 x 370 mm	2 (max. 4)
NBT203XL	180 l	560 x 665 x 945 mm	473 x 528 x 710 mm	3 (max. 6)
NBT203XL-28	280 l	700 x 690 x 1070 mm	603 x 580 x 790 mm	3 (max. 8)
NBT203XXL-1D	880 l	840 x 930 x 1890 mm	720 x 800 x 1530 mm	0 (max. 25)
NBT203XXL-2D	880 l	840 x 930 x 1890 mm	720 x 800 x 1530 mm	3 (max. 25)
NBT203XXL-5D	880 l	840 x 930 x 1890 mm	720 x 800 x 1530 mm	5 (max. 25)

► New Models

Model	Chamber Volume	External Dimensions (w x d x h)	Internal Dimensions (w x d x h)	Shelves
GNB-S203	42 l	408 x 482 x 550 mm	320 x 350 x 370 mm	2 (max. 4)
GNB-S203XL	180 l	560 x 665 x 945 mm	473 x 528 x 710 mm	3 (max. 8)
GNB-S203XXL	850 l	820 x 950 x 1840 mm	698 x 799 x 1528 mm	3 (max. 15)

NEW!

Shaking Incubators



Features

- Shakers work with a brushless DC motor for low noise level and reduced vibration
- Microprocessor controlled shaker and incubator
- Many optional accessories available such as platforms for flasks, tubes racks and 96-well plates or a universal sticky mat
- * Model with cooling function

CO₂ Incubators with Integrated Orbital Shakers

Model	Chamber Volume	External Dims. (mm) (w x d x h)	Internal Dims. (mm) (w x d x h)	Shelves	Shakers	Platform (cm)
NB203QS	179 l	560 x 665 x 945	473 x 528 x 710	1	1	30 x 33
NB206CL	179 l	560 x 665 x 945	473 x 528 x 710	-	2	30 x 33
NB206CX	450 l	820 x 780 x 1360	700 x 650 x 1030	-	2	52 x 52
NB206CXL	650 l	820 x 780 x 1740	700 x 650 x 1430	-	3	52 x 52
NB206CXXL	850 l	820 x 920 x 1840	700 x 800 x 1530	-	3	52 x 52

Incubators without CO₂ – for Benchtop, and Floor, also Stackable Units


Model	Chamber Volume	External Dims. (mm) (w x d x h)	Internal Dims. (mm) (w x d x h)	Shelves	Shakers	Platform (cm)
NB205	-	285 x 400 x 290	-	-	1	25 x 31
NB205Q	142 l	545 x 545 x 870	480 x 435 x 694	2	1	30 x 33
NB205V	-	600 x 660 x 995	-	-	1	45 x 45
NB205L	-	510 x 600 x 470	510 x 500 x 330	-	1	45 x 45
NB205LF *	-	550 x 790 x 555	510 x 500 x 330	-	1	45 x 45
NBT205	-	285 x 425 x 290	-	-	1	25 x 31
NBT205L	-	510 x 600 x 470	-	-	1	45 x 45
NBT205LF *	-	550 x 790 x 555	-	-	1	45 x 45
NBT207	149 l	600 x 870 x 730	510 x 505 x 577	1	1	45 x 45
NBT207F *	149 l	600 x 870 x 730	510 x 505 x 577	1	1	45 x 45

► New Models

Model	Chamber Volume	External Dims. (mm) (w x d x h)	Internal Dims. (mm) (w x d x h)	Shelves	Shakers	Platform (cm)
GNB-S206CL	180 l	820 x 780 x 1360	700 x 650 x 1030	-	2	52 x 52
GNB-S206CXL	650 l	820 x 775 x 1760	700 x 650 x 1430	-	3	52 x 52
GNB-S206CXXL	850 l	820 x 925 x 1860	700 x 800 x 1530	-	3	52 x 52
GNB-S205	-	285 x 400 x 290	-	-	1	25 x 31
GNB-T205P *	-	285 x 400 x 290	-	-	1	25 x 31
GNB-T205LP *	-	550 x 790 x 555	-	-	1	45 x 45
GNB-S207	150 l	600 x 870 x 730	-	1	1	46 x 46
GNB-T207P *	150 l	600 x 725 x 730	-	1	1	46 x 46

NEW!

Shakers



Features

- User-friendly and intuitive touch screen
- Programmable steps for temperature and speed, PID-controlled
- Units with separate remote controller can be used in incubators with a humidity up to 85 % and an ambient temperature of 4 °C to 60 °C, since the control box is separated from the shaker unit
- Brushless DC motor for low noise level and reduced vibration
- Shaking speed: 30 to 300 rpm (± 1 rpm); orbit diameter: 22 mm
- Timer: up to 99 h 59 min, or continuous operation
- NBT101S: optional platform to enable switch between orbital and reciprocal motion available
- NBT101MT: supplied with a basic platform with silicone rubber mat
- NBT101MC: consists of two platforms with mounted spring racks, which can be combined to one platform

► A sticky mat and various platforms for flasks, tube racks, 96-well microplates are available

Mini Shakers with Platform Size 30 x 33 cm

Cat. No.	Description
NB101S	Mini shaker
NBT101S	Mini shaker with touch screen
NB101SRC	Mini shaker with separate remote controller
NBT101SRC-V	Mini shaker with separate remote controller and vibration function
NBT101SRC	Mini shaker with separate remote controller and touch screen

Medium Shakers with Platform Size 46 x 45.5 cm

Cat. No.	Description
NB101M	Medium shaker
NBT101M	Medium shaker with touch screen
NBT101MRC	Medium shaker with separate remote controller and touch screen

Multi and Combination Shakers with Touch Screen and Platform Size 46 x 45.5 cm

Cat. No.	Description
NBT101MT	Multi shaker, orbital <u>or</u> reciprocal movement
NBT101MC	Combination shaker, orbital <u>and</u> reciprocal movement

Vacuum Evaporator

- Infrared glass lid spreads heat and infrared radiation to reduce time required for evaporation of samples using vacuum e.g. for DNA or RNA samples
- Temperature range: 35 °C to 65 °C (0.1 °C increments)
- 3 selectable rotation speeds: 1700, 2500, or 5000 rpm
- LED display for temperature and timer
- Vacuum pump and glass trap are built-in
- Easy control of sample without stopping spin by using the real scope
- Optional cold trap available



Cat. No.	Description
NB503CIR	Vacuum evaporator, dimensions (w x d x h): 213 x 255 x 225 mm

Water Bath



- Temperature range: 5 °C above ambient to 99 °C (± 0.5 °C at 37 °C)
- Convenient adjustment of water level with engraved level scale
- Heater located underneath the water chamber
- LED display for temperature, time, and power failure
- Supplied with flat bath cover
- Bath capacity: 10 l

Cat. No.	Description
NB301	Water bath, with flat cover, max. 10 l, dimensions (w x d x h): 270 x 345 x 225 mm
NB303	Orbital shaking water bath, without cover, max. 10 l, dimensions (w x d x h): 310 x 345 x 260 mm

Shaker Accessories

For each shaker, shaking incubator or incubator with integrated shaker, various support platforms and accessories are available. The catalogue number depends on the shaker model. Please contact us for further information.



2nd platform with 4 legs



platform with premounted spring wire rack



platform with premounted clamps



platform with premounted microplate rack

For more information please visit our website www.dunnlab.de
or contact us directly at info@dunnlab.de.