

Again available!

Terasaki Trays

- 60-well and 72-well trays
- Low profile, with 7.5 mm height
- For use in HLA typing, cytotoxicity, and other micro testing
- Sterile plates may be used for tissue culture
- Made of virgin medical-grade polystyrene
- Surface-treated for wettability
- Fit standard holders and can be used with manual and automated multiple syringe dispensers
- Approx. 10 µl volume per well
- Complete well-bottom is visible using 100X magnification
- Extreme flatness of plates minimizes the need to refocus between wells
- Plates are sterile
- Each tray is packaged with its own lid which has two tabs for easy orientation at opposite corners




Cat. No.	Description	Cs.
504S	60-well Terasaki plates, sterile	500
1004S	60-well Terasaki plates, sterile	1000
506S	72-well Terasaki plates, sterile	500
1006S	72-well Terasaki plates, sterile	1000

→ Please contact us for information on our other tissue typing products from Art Robbins Instruments.

Wieder im Programm!

Terasaki Platten

- 60-Well und 72-Well Terasaki Platten
- Niedrigprofilplatten mit einer Höhe von 7,5 mm
- Mögliche Anwendungen sind HLA, Zytotoxizitätsversuche und andere Mikrotest-Assays
- Sterile Platten können in der Zellkultur eingesetzt werden
- Aus hochwertigem, medizinischem Polystyrol
- Oberfläche der Platten ist feuchtigkeitsbehandelt
- Passen in Standardhalterungen und können mit manuellen und automatischen Mehrfachspritzendispensern verwendet werden
- Ca. 10 µl Volumen pro Well
- Der komplette Wellboden ist bei einer 100X-Vergrößerung immer sichtbar
- Durch große Ebenheit der Platten ist keine Neufokussierung nach jedem Well nötig
- Platten sind steril
- Jede Platte ist einzeln mit Deckel verpackt und hat zur einfachen Orientierung zwei Überhänge an gegenüberliegenden Enden



Kat. Nr.	Beschreibung	VE
504S	60-Well Terasaki Platten, steril	500
1004S	60-Well Terasaki Platten, steril	1000
506S	72-Well Terasaki Platten, steril	500
1006S	72-Well Terasaki Platten, steril	1000

→ Kontaktieren Sie uns für Informationen zu unseren anderen HLA-Produkten von Art Robbins Instruments.