

# Special N-Biotek Sale

## Up to 20 % Discount

### on Shakers, Incubators and Shaking Incubators

from May to July 2025

## Incubators



## Specifications

| Model        | Chamber Volume | External Dimensions (w x d x h) | Internal Dimensions (w x d x h) | Shelves     |
|--------------|----------------|---------------------------------|---------------------------------|-------------|
| NB203M       | 15,2 l         | 292 x 333 x 433 mm              | 224 x 200 x 340 mm              | 2 (max. 3)  |
| NB203        | 42 l           | 408 x 482 x 550 mm              | 320 x 350 x 370 mm              | 2 (max. 4)  |
| NB203XL      | 180 l          | 560 x 665 x 945 mm              | 473 x 528 x 710 mm              | 3 (max. 6)  |
| NB203XXL     | 850 l          | 820 x 950 x 1840 mm             | 698 x 799 x 1528 mm             | 3 (max. 15) |
| NBT203       | 42 l           | 408 x 482 x 550 mm              | 320 x 350 x 370 mm              | 2 (max. 4)  |
| NBT203XL     | 180 l          | 560 x 665 x 945 mm              | 473 x 528 x 710 mm              | 3 (max. 6)  |
| NBT203XL-28  | 280 l          | 700 x 690 x 1070 mm             | 603 x 580 x 790 mm              | 3 (max. 8)  |
| NBT203XXL-1D | 880 l          | 840 x 930 x 1890 mm             | 720 x 800 x 1530 mm             | 3 (max. 25) |

## Features

- Temperature range: 5 °C above ambient to 60 °C (± 0.25 °C at 37 °C)
- CO<sub>2</sub> range: 0 to 20 % (infrared sensor); CO<sub>2</sub> accuracy: ± 0.1 % at 5 %
- Stainless steel interior and perforated shelves
- Supplied with water tray for natural humidity
- Many optional accessories available such as additional shelves, access port, stacking kits, UV or hot air sterilization, Peltier cooling, oxygen control, and more!

## Standard Models

| Cat. No. | Description                                                       | List Price | Special Price |
|----------|-------------------------------------------------------------------|------------|---------------|
| NB203M   | MiniCell™ CO <sub>2</sub> portable incubator, 15.2 l              | 3,524 €    | 2,820 €       |
| NB203    | CO <sub>2</sub> incubator, 42 l                                   | 7,110 €    | 6,079 €       |
| NB203XL  | CO <sub>2</sub> incubator, 180 l                                  | 8,115 €    | 6,892 €       |
| NB203XXL | Large capacity CO <sub>2</sub> incubator with touch screen, 850 l | 19,702 €   | 16,272 €      |

## Models with Touch Screen

| Cat. No.     | Description                                                       | List Price | Special Price |
|--------------|-------------------------------------------------------------------|------------|---------------|
| NBT203       | CO <sub>2</sub> incubator with touch screen, 42 l                 | 7,980 €    | 6,783 €       |
| NBT203XL     | CO <sub>2</sub> incubator with touch screen, 180 l                | 8,738 €    | 7,396 €       |
| NBT203XL-28  | CO <sub>2</sub> incubator with touch screen, 280 l                | 9,932 €    | 8,363 €       |
| NBT203XXL-1D | Large capacity CO <sub>2</sub> incubator with touch screen, 880 l | 22,838 €   | 18,811 €      |


# Shaking Incubators



NBT207



NB203QS



NBT205



NB206CXXL



NBT207F

## Specifications

| Model     | Chamber Volume | External Dimensions (w x d x h) | Internal Dimensions (w x d x h) | Shelves | Shakers | Platform (cm) |
|-----------|----------------|---------------------------------|---------------------------------|---------|---------|---------------|
| NB203QS   | 179 l          | 560 x 665 x 945 mm              | 473 x 528 x 710 mm              | 1       | 1       | 30 x 33       |
| NB206CL   | 179 l          | 560 x 665 x 945 mm              | 473 x 528 x 710 mm              | -       | 2       | 30 x 33       |
| NB206CX   | 450 l          | 820 x 780 x 1360 mm             | 700 x 650 x 1030 mm             | -       | 2       | 52 x 52       |
| NB206CXL  | 650 l          | 820 x 780 x 1740 mm             | 700 x 650 x 1430 mm             | -       | 3       | 52 x 52       |
| NB206CXXL | 850 l          | 820 x 920 x 1840 mm             | 700 x 800 x 1530 mm             | -       | 3       | 52 x 52       |
| NB205     | -              | 285 x 400 x 290 mm              | -                               | -       | 1       | 25 x 31       |
| NBT205    | -              | 285 x 425 x 290 mm              | -                               | -       | 1       | 25 x 31       |
| NBT207    | 149 l          | 600 x 870 x 730 mm              | 510 x 505 x 577 mm              | 1       | 1       | 45 x 45       |
| NBT207F   | 149 l          | 600 x 870 x 730 mm              | 510 x 505 x 577 mm              | 1       | 1       | 45 x 45       |

## Features

- Temperature range: 5 °C above ambient to 60 °C (± 0.25 °C at 37 °C); NBT207F: 4 °C to 60 °C
- Shakers work with a brushless DC motor for low noise level and reduced vibration
- Timer: up to 48 hours, or continues operation
- Microprocessor controlled shaker and incubator
- Many optional accessories available such as platforms for flasks, tubes racks and 96-well plates or a universal sticky mat

## CO<sub>2</sub> Incubators with Integrated Orbital Shakers

| Cat. No.  | Description                                                        | List Price | Special Price |
|-----------|--------------------------------------------------------------------|------------|---------------|
| NB203QS   | CO <sub>2</sub> Incubator with 1 integrated orbital shaker, 179 l  | 10,862 €   | 9,116 €       |
| NB206CL   | CO <sub>2</sub> Incubator with 2 integrated orbital shakers, 179 l | 12,146 €   | 10,155 €      |
| NB206CX   | „aniCell™“, CO <sub>2</sub> incubator with 2 shakers, 450 l        | 24,801 €   | 20,400 €      |
| NB206CXL  | “aniCell”™, CO <sub>2</sub> incubator, 3 integrated shakers, 650 l | 28,345 €   | 23,269 €      |
| NB206CXXL | “aniCell”™, CO <sub>2</sub> incubator, 3 integrated shakers, 850 l | 33,246 €   | 27,327 €      |


## Benchtop Shaking Incubator and Stackable Shaking Incubator

| Cat. No. | Description                                                 | List Price | Special Price |
|----------|-------------------------------------------------------------|------------|---------------|
| NB205    | Mini shaking incubator                                      | 2,279 €    | 1,823 €       |
| NBT205   | Mini shaking incubator with touch screen                    | 2,476 €    | 1,981 €       |
| NBT207   | <b>NEW!</b> Stackable shaking incubator, 149 l              | 6,783 €    | 5,814 €       |
| NBT207F  | <b>NEW!</b> Stackable refrigerated shaking incubator, 149 l | 7,654 €    | 6,519 €       |


# Shakers



NBT101S inside an incubator



NB101S



NBT101M



NBT101MT



NBT101MC

## Features

- User-friendly and intuitive LCD touch screen
- Programmable steps for temperature and speed, PID-controlled
- Units with separate remote controller can be used in incubators with a humidity up to 85 % and an ambient temperature of 4 °C to 60 °C, since the control box is separated from the shaker unit
- Brushless DC motor for low noise level and reduced vibration
- Shaking speed: 30 to 300 rpm ( $\pm 1$  rpm); orbit diameter: 22 mm
- Timer: up to 99 h 59 min, or continuous operation
- NBT101S: optional platform to enable switch between orbital and reciprocal motion available
- NBT101MT: supplied with a basic platform with silicone rubber mat
- NBT101MC: consists of two platforms with mounted spring racks, which can be combined to one platform

► A sticky mat and various platforms for flasks, tube racks, 96-well microplates are available

## Mini Shakers with Platform Size 30 x 33 cm

| Cat. No.    | Description                                                        | List Price | Special Price |
|-------------|--------------------------------------------------------------------|------------|---------------|
| NB101S      | Mini shaker                                                        | 1,497 €    | 1,198 €       |
| NBT101S     | Mini shaker with touch screen                                      | 1,705 €    | 1,364 €       |
| NB101SRC    | Mini shaker with separate remote controller                        | 2,271 €    | 1,817 €       |
| NBT101SRC-V | Mini shaker with separate remote controller and vibration function | 3,044 €    | 2,436 €       |
| NBT101SRC   | Mini shaker with separate remote controller and touch screen       | 2,468 €    | 1,975 €       |

## Medium Shakers with Platform Size 46 x 45.5 cm

| Cat. No.  | Description                                                    | List Price | Special Price |
|-----------|----------------------------------------------------------------|------------|---------------|
| NB101M    | Medium shaker                                                  | 2,065 €    | 1,652 €       |
| NBT101M   | Medium shaker with touch screen                                | 2,221 €    | 1,777 €       |
| NBT101MRC | Medium shaker with separate remote controller and touch screen | 2,986 €    | 2,389 €       |

## Multi and Combination Shakers with Touch Screen and Platform Size 46 x 45.5 cm

| Cat. No. | Description                                                | List Price | Special Price |
|----------|------------------------------------------------------------|------------|---------------|
| NBT101MT | Multi shaker, orbital <u>or</u> reciprocal movement        | 2,608 €    | 2,087 €       |
| NBT101MC | Combination shaker, orbital <u>and</u> reciprocal movement | 3,144 €    | 2,516 €       |

## Vacuum Evaporator

- Infrared glass lid spreads heat and infrared radiation to reduce time required for evaporation of samples using vacuum e.g. for DNA or RNA samples
- Temperature range: 35 °C to 65 °C (0.1 °C increments)
- 3 selectable rotation speeds: 1700, 2500, or 5000 rpm
- LED display for temperature and timer
- Vacuum pump and glass trap are built-in
- Easy control of sample without stopping spin by using the real scope
- Optional cold trap available



| Cat. No. | Description                                                   | List Price | Special Price |
|----------|---------------------------------------------------------------|------------|---------------|
| NB503CIR | Vacuum evaporator, dimensions (w x d x h): 213 x 255 x 225 mm | 6,853 €    | 5,870 €       |

## Water Bath



- Temperature range: 5 °C above ambient to 99 °C ( $\pm 0.5$  °C at 37 °C)
- Convenient adjustment of water level with engraved level scale
- Heater located underneath the water chamber
- LED display for temperature, time, and power failure
- Supplied with flat bath cover
- Bath capacity: 10 l


| Cat. No. | Description                                                                        | List Price | Special Price |
|----------|------------------------------------------------------------------------------------|------------|---------------|
| NB301    | Water bath, with flat cover, max. 10 l, dimensions (w x d x h): 270 x 345 x 225 mm | 957 €      | 766 €         |

## Shaker Accessories

For each shaker, shaking incubator or incubator with integrated shaker, various support platforms and accessories are available. The catalogue number depends on the shaker model. Please contact us for further information.



2<sup>nd</sup> platform with 4 legs



platform with premounted spring wire rack



platform with premounted clamps



platform with premounted microplate rack

For more information please visit our website [www.dunnlab.de](http://www.dunnlab.de)  
or contact us directly at [info@dunnlab.de](mailto:info@dunnlab.de).