



MADE IN 
USA



ROLLER APPARATUS

Keep your research rolling

One size fits most.
CUSTOM fits the rest.



*Most Bellco products can be configured to suit your unique needs.



ROLL WITH PRECISION. GROW WITH CONFIDENCE.

Bellco's Roller Apparatus fits seamlessly in incubators and scales with your workflow. With customizable spacing, non-slip chemically resistant rollers, digital touchscreen control, and connectivity via Wi-Fi or Ethernet, it's precision-made for high-yield cell culture.

DESIGNED FOR HOW YOU WORK

Roller Apparatus can be designed and engineered to fit your needs and requirements. Just ask us!



MODELS
Modular (easy to expand)
Fixed Frame

NUMBER OF BOTTLES
1 to 176
Custom Quantities Available

ROLLER BOTTLE VOLUMES
670, 850, 1170, 1330, & 1585mL

COMMUNICATION
Ethernet
Wi-Fi
4-20mA
RS-485

LOCATION
Bench Top
Floor Mount
with Casters

VOLTAGE
100V
115V
230V

CONTROLLER
Digital Display
PID with Digital Display
PID with Touchscreen Interface

OPTIONS
Battery Back-up
Temperature Display
Motion Alarm
Data Logger

NUMBER OF DECKS
1 to 11 position

DECK SPACING
6.375in
7.125in
Custom Spacing

ROLLERS
Chemically Resistant
Non-Slip EDPM

RPM
0.1 to 4.0rpm
3.0 to 9.0rpm
Custom Settings

DRIVES LOCATIONS
Top Mount
Bottom Mount
Side Mount

ROLLER APPARATUS

Made for each other. Bellco Incubators and Roller Apparatus.



Bellco's Roll-In Incubators are built to accommodate large-scale roller apparatus setups with ease. Designed for uniform temperature distribution and optimal air-flow, these units maintain a consistent environment for high-yield cell culture. With durable stainless-steel interiors, easy-access doors, and reliable controls, Bellco incubators make scaling up as simple as rolling in.

DESIGNED FOR HOW YOU WORK

Most Bellco products
can be configured to
suit your unique needs.



CONTACT US:

Tel.: +49 (0) 26 83 / 4 30 94

E-mail: info@dunnlab.de

Internet: www.dunnlab.de

